



Lab #2 1/23/14

1. Set up rotating stage, photodetector, & optical path



2. Plugged stage in & connected PD to scope.
3. Placed polarizer in path to determine polarization of beam.
4. Moved stage off track for better angle sweeping & realigned



5. Placed iris in path to make beam smaller for PD
6. Placed polarizer between  $M_1$  &  $M_2$
7. Polarizer at  $261^\circ$  - almost no reflection (polarizer)  
No reflection again at  $81^\circ$

8. Light intensity at glass angle  $56^\circ$



800 mV average signal from PD (background  $\sim 800$  mV)

Lights off:

Polarizer @  $261^\circ$  (perfect extinction)

Background: 0V

<u><math>\theta_{\text{glass}}</math></u>	<u>PD Voltage</u>
$56^\circ$	0 mV
$54^\circ$	150 mV
$52^\circ$	250 mV
$50^\circ$	650 mV
$45^\circ$	1600 mV
$40^\circ$	2800 mV
$35^\circ$	4000 mV
$30^\circ$	5000 mV
$25^\circ$	5200 mV
$20^\circ$	5000 mV
$15^\circ$	5200 mV
$10^\circ$	5400 mV

~~5800 mV~~

Experiment #2

Next we did the other polarization

$$\phi_{\text{polarizer}} = 171^\circ$$

<u><math>\theta_{\text{glan}}</math></u>	<u>PD Voltage</u>
56°	11,300 mV
54°	11,300 mV
52°	

Realized PD was at  $V_{\text{sat}}$ , placed resistor (brown black black)

<u><math>\theta_{\text{glan}}</math></u>	<u>PD Voltage</u>
56°	200 mV
54°	150 mV
52°	150 mV
50°	50 mV

Needed ~~different~~ ~~smaller~~ resistor:  $10^3 \text{ k}\Omega$

<u><math>\theta_{\text{glan}}</math></u>	<u>PD Voltage</u>
56°	330 mV
54°	300 mV
52°	260 mV
50°	240 mV
45°	200 mV
40°	165 mV
36°	137 mV
30°	120 mV
25°	105 mV
20°	93 mV
15°	81 mV
10°	77 mV


## Experiment #2

$$\theta_p = 261^\circ \leftarrow \text{natural}$$

$$R = 10^3 \text{ k}\Omega$$

Photodetector  
on varifocall  
glass for  
transmission

$\theta_{\text{glass}}$	PD Voltage
$56^\circ$	2600 mV
<del>540</del>	2660 mV
$52^\circ$	2740 mV
$50^\circ$	2740 mV
$45^\circ$	2700 mV
$40^\circ$	2740 mV
$35^\circ$	2580 mV
$30^\circ$	2500 mV
$25^\circ$	2500 mV
$20^\circ$	2450 mV
$15^\circ$	2440 mV



$$\theta_p = 171^\circ \text{ max}$$

$\theta_{\text{glass}}$	PD Voltage
$56^\circ$	1850 mV
$54^\circ$	2100 mV
$52^\circ$	2100 mV
$50^\circ$	2180 mV
$45^\circ$	2200 mV
$40^\circ$	2420 mV
$35^\circ$	2500 mV
$30^\circ$	2500 mV
$25^\circ$	2500 mV
$20^\circ$	2500 mV
$15^\circ$	2580 mV

### Experiment #3

Plastic slab thickness: 1.7 cm

Deflection: 1.3 cm

### Experiment #4

~~Exp.~~  $f_1 = 24.4 \text{ mm}$  } lenses  
 $f_2 = 200 \text{ mm}$

Furthest: 1.75 m

Next: 49 cm